




**WICKLAND MEDIKA IMANTARA**  
MEDICAL EQUIPMENT , AESTHETICS AND LABORATORY SUPPLIER

**UROLOGY CATALOG**



 **LOCATION**  
PALEM SEMEMI RAYA NO. 12-12A  
KEL. SEMEMI, KEC. BENOWO, SURABAYA  
JAWA TIMUR, INDONESIA

 **CONTACT US**  
031-33300312  
031-33100412

 **EMAIL**  
WICKLAND MEDIKA@GMAIL.COM

Single Use Grasping Forceps



● Indication :

The disposable grasping forceps is indicated to grasp stents, stone or other foreign body objects during flexible cystoscopy.



Single Use Grasping Forceps

● Features :



Large jaw design to take larger scale of foreign body.



The head is precisely designed for grasping foreign body accurately and firmly.




The end of the jaw is smooth after elaborately grinding, without damage to the endoscopy.

● Product Model :

| Model | Order Number  | Outer Diameter (mm) | Working Length (mm) | Working Channel (mm) | Coating  | Jaw Type                |
|-------|---------------|---------------------|---------------------|----------------------|----------|-------------------------|
| C     | MD-C-GF1806 ★ | 1.8                 | 600                 | 2.0                  | Optional | Rat tooth Alligator jaw |
|       | MD-C-GF1812   |                     | 1200                |                      |          |                         |

# Single Use Urinary Guidewire



● Indication :

Used to gain ureteral access to establish a tract and to assist during urological procedures.



# Single Use Urinary Guidewire


● Features :

- ⌘ Super smooth hydrophilic surface ensures excellent control and direction selection without any damage to tissues.
- ⌘ Yellow and black coating provides high visualization under X-ray.
- ⌘ Anti-twist Nitinol core wire is easier to pass through complicated working channel.
- ⌘ Conical tip design makes precise control and positioning.
- ⌘ Superfine spring inside with great toughness ensures easy access.

● Product Model :

| Model                               | Order number   | Outer diameter (inch) | Working length (cm) | Tip      |
|-------------------------------------|----------------|-----------------------|---------------------|----------|
| A<br>(Hydrophilic Tip & Zebra Body) | MD-A2-UW3515 ★ | 0.035inch (0.89mm)    | 150                 | Straight |
| B<br>(Hydrophilic)                  | MD-B2-UW3215 ★ | 0.032inch (0.81mm)    |                     | Straight |

Single Use Ureteral Access Sheath



● Indication:

Used as a working channel to maintain a previously established access, and to assist the flexible endoscope and other equipment into the urinary tract.



Single Use Ureteral Access Sheath

● Features:

- ✘ Hydrophilic coated sheath creates a low-friction surface as the ureter is dilated.
- ✘ Strengthening wire makes endoscope and instruments enter smoothly, and the stainless steel wire layer provides better compression resistance.
- ✘ The tapered tip can be gradually expanded into place smoothly.
- ✘ The surface scale easily control the entry distance.

● Product Model :

| Model | Order Number | I.D.(Fr) | O.D.(Fr) | Working Length(mm) |
|-------|--------------|----------|----------|--------------------|
| A     | MD-A-AS-1045 | 10       | 12       | 450                |
|       | MD-A-AS-1055 | 10       | 12       | 550                |
|       | MD-A-AS-1245 | 12       | 14       | 450                |
|       | MD-A-AS-1255 | 12       | 14       | 550                |
| B     | MD-B-AS-1035 | 10       | 12       | 350                |
|       | MD-B-AS-1045 | 10       | 12       | 450                |
|       | MD-B-AS-1235 | 12       | 14       | 350                |
|       | MD-A-AS-1245 | 12       | 14       | 450                |

Disposable Ureteral Stent



● Indication:

Used for temporary internal drainage from the ureteropelvic junction to the bladder.



Disposable Ureteral Stent

● Features:

- ✘ Temperature-controlled disposable ureteral stent, it becomes soft into the body for 3-5 minutes to reduce mucosal stimulation;
- ✘ Resilient, smooth, it can be kept in the body for one month;
- ✘ Medical high polymer materials with superb biocompatibility;
- ✘ Big lumen with multiple holes to improve smooth flow;
- ✘ one side open , both sides pen, both sides close for your selection;


● Product Model:

| Model                     | Order number | Outer diameter (Fr) | Working length (mm) |
|---------------------------|--------------|---------------------|---------------------|
| C<br>Single open Double J | MD-C-JC-5F   | 5                   | 260                 |
|                           | MD-C-JC-6F   | 6                   |                     |
| D<br>Double open double J | MD-D-JC-4F   | 4                   |                     |
|                           | MD-D-JC-5F ★ | 5                   |                     |
|                           | MD-D-JC-6F ★ | 6                   |                     |
|                           | MD-D-JC-7F   | 7                   |                     |
|                           | MD-D-JC-8F   | 8                   |                     |

Remark: stent package is optional with pusher, clamp, PTFE guidewire.

# Single Use Urinary Stone Retrieval Basket

# Single Use Urinary Stone Retrieval Basket



● **Indication:**

Used in an endoscopic procedure for removal of urinary stones or foreign objects from the ureter or the kidney.

● **Features:**

- ⌘ Advanced Nitinol wires ensures the excellent flexibility and memory, keeping the shape after extraction.
- ⌘ Multi-model customization meets the requirements of doctors.
- ⌘ Unique catheter design, providing excellent propulsion and hand feeling.
- ⌘ Sterile & single use.

● **Product Model:**

| Model | Order number      | Maximum width (Fr) | Basket width (mm) | Working length (cm) | Shape           |
|-------|-------------------|--------------------|-------------------|---------------------|-----------------|
| C     | MD-C-UN05908115   | 1.7Fr (0.59mm)     | 8                 | 115                 | 3 Wires Tipless |
|       | MD-C-UN07311115 ★ | 2.2Fr (0.73mm)     | 11                |                     | 3 Wires Tipless |
| E     | MD-E-UN05908115   | 1.7Fr (0.59mm)     | 8                 |                     | 4 Wires Tipped  |
|       | MD-E-UN07310115   | 2.2Fr (0.73mm)     | 10                |                     | 4 Wires Tipped  |
| F     | MD-F-UN05908115   | 1.7Fr (0.59mm)     | 8                 |                     | 4 Wires Tipless |
|       | MD-F-UN07310115 ★ | 2.2Fr (0.73mm)     | 10                |                     | 4 Wires Tipless |



# Single Use Urological Retrieval Coil



● **Indication :**

Used to prevent stone migration to the kidney and to facilitate extraction of stone fragments.



# Single Use Urological Retrieval Coil


● **Features :**

- ✕ The guiding tip is flexible and can pass through stones smoothly.
- ✕ The retrieval coil is made of nitinol core wire with strong elasticity, which can seal the stone tightly, effectively reduce the displacement of the stone, and improve the stone clearance rate.
- ✕ The surface is covered with hydrophilic coating, which is easy to pass, reducing ureteral damage.
- ✕ The external control handle is convenient for procedure.

● **Product Model :**

| Model | Order number   | Sheath O.D.<br>(Fr) | Coil O.D.<br>(mm) | Working length<br>(cm) |
|-------|----------------|---------------------|-------------------|------------------------|
| B     | MD-B-RC-1012 ★ | 3                   | 10                | 120                    |
|       | MD-B-RC-1010   | 3                   | 10                | 100                    |

Single Use Biopsy Forceps



● Indication:

This device is used to enter the gastrointestinal tract via endoscope to obtain tissue samples for pathology.



Single Use Biopsy Forceps


● Features:

- ✕ Large cup design to get more tissues.
- ✕ Rigid four-link structure to make the sampling easier.
- ✕ Sharp cutting to get the right amount of tissue, with less bleeding.
- ✕ Sterilized package to prevent the risk of cross-infection.

● Product Model:

| Model | Order Number  | Outer Diameter (mm) | Working Length (mm) | Working Channel (mm) | Needle | Jaw Type |
|-------|---------------|---------------------|---------------------|----------------------|--------|----------|
| A     | MD-A-BF1106   | 1.0                 | 600                 | 1.2                  | No     | Oval     |
|       | MD-A-BF1806   | 1.8                 | 600                 | 2.0                  |        | Oval     |
|       | MD-A-BF1112   | 1.0                 | 1200                | 1.2                  |        | Oval     |
|       | MD-A-BF1812   | 1.8                 | 1200                | 2.0                  |        | Oval     |
| C     | MD-C-BF1106   | 1.0                 | 600                 | 1.2                  |        | Serrated |
|       | MD-C-BF1806   | 1.8                 | 600                 | 2.0                  |        | Serrated |
|       | MD-C-BF1112 ★ | 1.0                 | 1200                | 1.2                  |        | Serrated |
|       | MD-C-BF1812 ★ | 1.8                 | 1200                | 2.0                  |        | Serrated |

Single Use Urethral Dilator



● Indication :

Used for dilation of the ureter prior to uretero renoscopy and also for sequential dilation of the intramural tunnel.



Single Use Urethral Dilator

● Features :

- ✕ Hydrophilic coating can effectively reduces ureteral damage.
- ✕ Uniform taper to reduce trauma.
- ✕ Biocompatible material has good hardness and toughnes,improving comfort.

● Product Model :

| Model                              | Order number   | O.D. (Fr)                                | Working length (mm) |
|------------------------------------|----------------|--|---------------------|
| A kit<br>Straight shape for female | MD-A-UD-0060 ★ | 6, 8, 10, 12, 14, 16, 18, 20, 22, 24     | 600                 |
| B kit<br>S shape for male          | MD-B-UD-0040 ★ | 6, 8, 10, 12, 14, 16, 18, 20, 22, 24, 26 | 400                 |



**WICKLAND MEDIKA IMANTARA**  
MEDICAL EQUIPMENT , AESTHETICS AND LABORATORY SUPPLIER

**CONTACT US:**



031-33300312  
031-33100412



WICKLANDMEDIKA@GMAIL.COM