

# Endocam FHD All in One WKL - 21



## SPEISIFICATION

Camera head	<ul style="list-style-type: none"> <li>Full HD imaging (1920 x 1080)</li> <li>Resolution CMOS1,000,000 Pixels high resolution</li> </ul>
Frequency	<ul style="list-style-type: none"> <li>50/60 Hz</li> </ul>
Device protection	<ul style="list-style-type: none"> <li>Emi Filter, 2x Fuse, 3A / 250V</li> </ul>
Video system	<ul style="list-style-type: none"> <li>Image resolution: FHD 1920 x 1080</li> <li>Scanning system: Progressive Scanning</li> <li>White balance automatically</li> </ul>
Recording (optional)	<ul style="list-style-type: none"> <li>Image format: JPG 1920 x 1080</li> <li>Video format: Full HD 1920 x 1080</li> <li>MPEG-4 (H.264 Codec); Image rate 60fps</li> <li>Connection: USB</li> </ul>
Connections	<ul style="list-style-type: none"> <li>Video out terminal</li> <li>1x HDMI</li> </ul>
Protection class	<ul style="list-style-type: none"> <li>I</li> </ul>
Safety class	<ul style="list-style-type: none"> <li>BF</li> </ul>
Standard	<ul style="list-style-type: none"> <li>SNI IEC 60601-1:2014</li> <li>IEC 60601-2-18:2009</li> </ul>
Features	<ul style="list-style-type: none"> <li>White balance</li> <li>Freeze function</li> <li>Brightness, contrast, saturation adjustment</li> <li>Zoom in / Zoom out</li> <li>Touch Buttons</li> <li>USB storage</li> <li>Support connects same series</li> <li>Photo capture / Video recording</li> </ul>
Language	<ul style="list-style-type: none"> <li>Indonesian/Chinese/English/Japanese/Korean/Spanish/German/Russian</li> </ul>
Weight	<ul style="list-style-type: none"> <li>12,5 Kg</li> </ul>