



WICKLAND MEDIKA IMANTARA
MEDICAL EQUIPMENT , AESTHETICS AND LABORATORY SUPPLIER



LOCATION

PALEM SEMEMI RAYA NO. 12-12A
KEL. SEMEMI, KEC. BENOWO, SURABAYA
JAWA TIMUR, INDONESIA




CONTACT US
031-33300312
031-33100412



EMAIL
WICKLANDMEDIKA@GMAIL.COM

KL-BC102

Electric Hospital Bed
5-Functions



Technical Parameters


External (L*W*H):	2240*1030*450-700mm
Safe load:	240kg
Back-rest Adjustment:	0-75°
Knee-rest Adjustment:	0-35°
Reverse Trendelenburg:	0-12°
Trendelenburg:	0-12°

Technical Configuration

Silent castors:	4 pcs
Control box:	1 set
Guardrail:	1 set
ABS drainage hooks:	2 pcs

KL-BC103

Electric Hospital Bed
5-Functions



Technical Parameters


External (L*W*H):	2150*1065*(540-830)mm
Mattress Platform:	1925*900mm
Safe Working Load:	240kg
Back-rest Adjustment:	0-75°(±10°)
Knee-rest Adjustment:	0-35°(±10°)
Reverse Trendelenburg:	0-12°
Trendelenburg:	0-12°

Technical Configuration

Silent castors:	4 pcs
Imported motor:	1 set
Guardrail:	1 set
Central locking pedal:	1 set

KL-BC104

Electric Hospital Bed
3-Functions



Technical Parameters


External (L*W*H):	2150*1035*450-730mm(±10mm)
Input power:	120VA
Safe Working Load:	240kg
Back-rest Adjustment:	0-75°
Knee-rest Adjustment:	0-35°
Height Adjustment:	280mm

Technical Configuration

Silent castors:	4 pcs
Push tod motor :	3 pcs
Guardrail :	1 set
Head and foot board :	1 set

KL-BC105

Electric Hospital Bed
2-Functions



Technical Parameters


Length:	2150mm
Width:	960mm
Height:	500(±10mm)
Mattress Platform:	1930*900mm
Safe Working Load:	240kg
Back-rest Adjustment:	0-75°
Knee-rest Adjustment:	0-35°

Technical Configuration

Manual crank:	2 pcs
Silent castors:	4 pcs
Guardrail:	1 set

KL-BC108

Manual Hospital Bed
3-Functions



Technical Parameters

External (L*W*H):	2150*960*450-730(±10mm)
Mattress Platform:	1930*900mm
Safe Working Load:	240kg
Back-rest Adjustment:	0-75°
Knee-rest Adjustment:	0-35°
Material:	stainless steel/ABS

Technical Configuration

Manual crank:	3 pcs
Silent castors:	4 pcs
Guardrail:	1 set

KL-BC109

Manual Hospital Bed
5-Functions



Technical Parameters


External (L*W*H):	2150*980*450-720mm
Mattress Platform:	1930*900mm
Safe Working Load:	240kg
Back-rest Adjustment:	0-75°
Knee-rest Adjustment:	0-35°
Reverse Trendelenburg:	0-12°
Trendelenburg:	0-12°

Technical Configuration

Manual crank:	4 pcs
Silent castors:	4 pcs
Guardrail:	1 set

KL-BC110

Electric Gynecological Birth Bed
5-Functions



Technical Parameters


Bed (L*W*H):	2050mm*980mm*680-900mm
Working environment temperature:	10 ° -40 °
Total input power:	400VA
Safe working load:	170KG
Back panel folding angle:	70 °
Overall upward tilt angle:	17 °
Overall downward tilt angle:	17 °

Technical Configuration

Handle controller:	1 pc
Leg support:	2 pcs
Movable panels:	1 pc
Guardrail:	1 set
Silent castors:	4 pc
Stainless steel basin:	1 set

KL-BC111

Manual Orthopedics Traction Bed
5-Functions



Technical Parameters

External (L*W*H):	2160*1060*460-720mm
Mattress Platform:	1925*900mm
Safe Working Load:	250kg
Back-rest Adjustment:	0-75°(±10°)
Knee-rest Adjustment:	0-35°(±10°)
Reverse Trendelenburg:	0-12°
Trendelenburg:	0-12°

Technical Configuration

Guardrail:	1 set
Head and foot board:	1 set
Silent castors:	4 pcs
Stainless steel drawing frame:	1

KL-BC129

Flat bed

Technical Parameters

Product size(L*W*H):	2000*960*500(±10mm)
Bed size:	1930*900mm
Headboard:	can be customized



KL-BC130

Flat bed

Technical Parameters

Product size(L*W*H):	2000*960*500(±10mm)
Bed size:	1930*900mm



KL-BC131

Single rocking nursing bed

Technical Parameters

Product size(L*W*H):	2000*960*500(±10mm)
Bed size:	1930*900mm
Back Angle:	0-75°
Tabula:	can be customized
Guard rail:	can be customized



KL-P001

Mattress



Technical Configuration & Optional

Mattress size: 190mm*90mm
6cm ☐ 8cm ☐ 10cm ☐ 12cm ☐ 14cm ☐
Mattress size & thickness: can be customized

KL-P002

Mattress



Technical Configuration & Optional

Mattress size: 190mm*90mm
6cm ☐ 8cm ☐ 10cm ☐ 12cm ☐ 14cm ☐
Mattress size & thickness: can be customized

KL-P003

Air Mattress



Technical Configuration & Optional

Mattress size: 190mm*90mm
Applicable to a variety of hospital scenarios

KL-P004

AirMattress



Technical Configuration & Optional

Mattress size: 190mm*90mm
Applicable to a variety of hospital scenarios

KL-LJ045

Hanger



KL-LJ045-1

IV Pole



Optional:

Material: Stainless steel
Outer tube thickness: 19mm
Inner tube thickness: 16mm


KL-HL001

Guard Rail



KL-HL002

Rail



KL-CA011

Castor



KL-CA012

Belt Of Equipment



KL-HL003

Headboard



KL-HL004

Headboard



Suitable for all kinds of beds,
please consult for more information



*Blowing Procees
Suitable for all kinds of beds,
please consult for more information



KL-CG081

Bedside Table

Size: 480*465*745mm (L*W*H)



KL-CG082

Bedside Table

Size: 480*480*760mm(L*W*H)



Option:



KL-CG083

Bedside Table

Size: 480*465*785mm (L*W*H)



KL-CG084

Bedside Table


Size: 480*400*790mm(L*W*H)



KL-CG080

Bedside Table

Size: 480*460*900mm(L*W*H)



KL-CG085

Bedside Table

Size: 480*480*760mm(L*W*H)



KL-CG086

Bedside Table

Size: 480*480*760mm(L*W*H)



KL-CG087

Bedside Table

Size: 480*480*760mm(L*W*H)



KL-CZ301

Overbed Table
Size:940*450*760-1040mm
For more information, please contact us.



KL-CZ302

Overbed Table
Size:940*450*760-1040mm



KL-CZ303

Overbed Table
Size:765*385*700-1010mm



KL-CZ304

Dining table
Size:870-1020*350*28mm




KL-CZ305

Dining table
Size:1030*300*28mm



KL-CL050

Curtain&Track



Track styles:	U-shaped, straight, L-shaped, etc.
Track material:	Aluminum alloy
Material:	Polyester
Height:	2.8m
Weight:	200g/m²-250g/ m²
Color & Track size:	customizable

KL-CL051

Fold Screen



Single screen size:	50*36*180cm
Frame material:	stainless steel
Material:	fiber
Mute universal castors:	1 set
Detachable shade cloth:	1 set
Stainless steel frame:	1 set

KL-LJ046

IV Pole

Product size: 025*39.5cm
(more sizes are optional)



KL-LJ047

Trash Can

Product size:
25*39.6cm(12L)
29*43.5cm(20L)



KL-LJ048

Kick Bucket

Product size:
33.5*48cm





KL-EA02

Examination Lamp

Technical Configuration & Optional

Material:	Stainless steel	Light Angle:	24°
The diameter:	56cm	Working voltage:	AC180-260V
Wheel:	6cm	Working current:	600MA/W ±5%
Illumination:	1750LM(30000Lux)	Power:	1*20W



KL-WS001

Weighing scale

Technical Configuration & Optional

Dimension:	62x29.5x106cm
N.W.:	11.33kg
G.W.:	13.4 kg
Packing size:	107x36x28cm(1pcs/ctn)
Use 4XAA Batteries:	(Not included)



KL-WS002

Weighing scale

Technical Configuration & Optional


Dimension:	62x29.5x106cm	Use 4XAA Batteries
N.W.:	11.33kg	(Not included)
G.W.:	13.4 kg	
Packing Volume:	107x36x28cm	
For more details, please consult		

KL-BD101

Vision chart lightbox

Technical Configuration & Optional

Product size:	619*304*25mm
Light source:	SMD2835 LED96
Color temperature:	6000k
Voltage:	AC90V-240V
Frequency:	50/60Hz
Installation method:	Wall-mounted
Observation screen brightness:	≥90%



X-LED Series

X ray film viewer

Technical Configuration & Optional

Model	X-Ray viewerfilm 1 plage
Frame size(L*W*H):	472X560X17mm
Visual size(L*W*H):	366X455
Energyconsumption:	30 W
Color temperature:	8000K
Input voltage:	Integrated power supply100-240V.50/60Hz
luminance:	Standard: high 4000cd/cm'(3000Lux)



KL-FS002

Footstool

Technical Configuration & Optional

Material:	Stainless steel
Product size:	47*22.5*23cm
More options please consult	




KL-ZY001

Treat-waiting Chair



Material: PU/PVC

Material: stainless steel



Technical Configuration & Optional

Product size:	180*74*80cm
Number of seats:	3 (customizable)
PU cushion:	304 stainless steel
Seat thickness:	1.5cm
Color:	optional
Main material:	customizable

KL-ZY002

Nurse Stool



Technical Configuration & Optional

Product size:	35*8*51-67cm
Cushion thickness:	8cm
Cushion diameter:	35cm
Cushion adjustable height:	51-67cm
Product weight:	4.8kg
Color:	Customizable

KL-ZY003

Office Chair



Technical Configuration & Optional

Product size:	55*56*103cm
Sitting height:	45cm
Backrest:	1
Armrest:	2 pcs
Bow leg:	1 set
Seat:	1

KL-ZY004

Attendant chair



Technical Configuration & Optional

Seat size (L*W*H):	640*720*940 mm	2-inch castors:	4
Bed size(L*W*H):	640*1900*420 mm	Non-slip feet:	4
Main material:	carbon steel	Pillow:	1
Arm:	1 set	Color:	customized

KL-ZY005

Transfusion Chair



Technical Configuration & Optional

Seat Size(L*W*H):	1150*800*1200mm	IV Pole:	1 pc
Backrest to ground height:	1200mm	Utility Basket:	1 pc
Cushion thickness:	60mm	Sponge size:	7cm
IV Pole Previsions:	2pcs	Color:	customized

KL-ZY006

Luxury Attendant Chair



Technical Configuration & Optional

Seat size:	710*720*950mm
Bed size:	710*1940*420mm
Material:	PU/wood
Wood Arm Frame:	1 set
2" Castors:	2 pcs
Anti-slip Feet Cover:	2 pcs
Pillow:	1 pc

KL-ZY007

Electric Dialysis Chair

Technical Configuration & Optional

Size(chair):	1000×600×1300mm
Size(open):	1880×880×600mm
Carton Size:	970x820x750mm
N.W:	44kg
G.W:	50kg
Seat cushion thickness:	120mm
Control handle:	1
Castor:	4
Color:	customizable



KL-ZY008

Escort chair

Technical Configuration & Optional

Size(chair):	1100x770x1280mm
Carton:	1pcs
Table size :	480*260mm
Carton Size:	970x820x750mm
N.W:	44kg
G.W:	50kg
Weight limit:	120KG
5" castor:	4
Color:	customizable



KL-ZY009

Companion chair

Technical Configuration & Optional

Product size:	830*780*930mm
Unfold size:	830*2000*930mm
PUcushion:	PU
Posture:	Sit / lie down
Number of seats:	1 (customizable)
Color:	customizable



KL-BC114

Child's cot



Technical Parameter

External(L*W*H):	1880*900*1260mm
Bed Height:	585mm
Mattress Platform:	1770*720mm
Back-rest Adjustment:	0-75°(±10°)
Knee-rest Adjustment:	0-35°(±10°)

Optional

Castors:	4 pcs
Color:	Blue ☒ Pink ☐

Technical Configuration

Manual Cranks:	2pcs
5" Covered Castors:	4pcs
Full Length Side Rail:	1set
Plastic Bed Ends with Cartoon Bear:	1pc
IV Pole Prevision:	4pcs
IV Pole:	1pc
Drainage Hook:	2pcs
BufferAir Spring:	4pcs

KL-BC115

Child's cot



Technical Configuration & Optional

External size(L*W*H):	800*500*800(mm)
Material:	Stainless steel
Color:	can be customized

KL-BC116

Baby's crib



Technical Configuration & Optional

External size(L*W*H):	810*520*990 mm
Trendelenburg:	18°
Bassinet:	1pc
3" Covered Castors:	4pcs
IV Pole:	1pc
Adjustable Air Spring:	4pcs
Color:	can be customized

KL-BC117

Baby's crib



Technical Configuration & Optional

Trendelenburg:	0-12°(±5°)
ReversedTrendelenbury:	0-12°(±5°)
Bassinet:	1pc
3"covered Castors:	4pcs
Size:	875*525*(650-900)mm
Mattress size:	2cm
Load Capacity:	>30kg
Packing Size:	910*550*850 mm
	1 pc/carton
Color:	can be customized



Animation Demonstration

KL-2F

Electric Gynecological Bed



Technical Parameters

Produce size (L*W):	1950±50mm*700±20mm
Bed surface:	(650-920mm)±50mm
The table height adjustment:	270±50mm
Forward tilt:	≥25°
Backward tilt:	≥25°
Fold on the back panel:	≥70°
Fold down on the back panel:	≥5°
Outward folding ofleg plate:	≥90°
Input power:	150VA±10%
Power supply voltage:	AC 220V±22V 50Hz±1Hz

Leg support:	1 pc	Foot remote control:	1 set
Kick board:	1 pc	Handheld remote control:	1 set
Handle rod:	1 pc	Stainless steel dirt basin:	1 set
Power cord:	1 set	Color:	Blue ☐ Green ☐ Red ☐ Pink ☐
Arm rest:	1 pc		

KL-FA.IV

Electric Gynecological Bed



Technical Parameters

Product size:	140*58*(97-56)cm
Forward tilt:	≤70°
Backward tilt:	≤16°
Angle between back and seat:	≤43°
Auxiliary plate:	optional
Castor:	4 pcs
Color:	Blue ☐ Green ☐ Red ☐ Pink ☐

KL-2A

Gynecological Examination Table



Technical Parameters

Product size:	1850*650*820±20mm
Back board upwards:	50°±5°
Leg board downwards:	0-90°



Optional

Color:	customizable
--------	--------------

KL-BC122

Gynecological Examination Table



Technical Parameters

External (L*W*H):	1880*600*680mm
Back-rest Adjustment:	0-70°
Drawers size :	540*400*100mm
Sliding stool size:	400*320*120mm
Stainless steel basin:	500*350*48mm
Packing size:	1440*740*770mm
Safe Working Load:	250kg




Technical Configuration & Optional

Mattress:	1 pc	Color of mattress:	customizable
Leg supports:	2 pcs		
Stainless steel drain basin:	1 pc		
Stool:	1 pc		

KL-3A

Gynecological Examination Table



Technical Parameters

Product size:	1850*650*820+20mm
Back board upwards:	0-85°
Leg board downwards:	0~90°
Carrying capacity:	150kg

Optional

Color:	customizable
--------	--------------

KL-BC123

Examination Table



Technical Parameters

External (L*W*H):	1800*600*600mm
Bed body:	1 set

Optional

Castors:	optional
Color:	customizable

KL-BC121

Examination Table



Technical Parameters


External (L*W*H):	1900*600*680mm
Back-rest Adjustment:	0-60°(±10°)
Bed body:	1

Optional

Castors:	optional
Color:	customizable

KL-BC120

Examination Table



Technical Parameters

External Length:	1940mm
External Width:	700mm
Table Height Min.:	490mm
Table Height Max.:	890mm




Optional

Castors:	optional
Color:	customizable

KL-BC127

Examination Table



Technical Parameters

External Length:	1700mm
External Width:	600mm
Table Height Min.:	690mm
Table Height Max.:	980mm



Optional

Castors:	optional
Color:	customizable

KL-BC125

Examination Table



Technical Parameters

External (L*W*H):	1900*650*700mm
Back-rest Adjustment:	60°
Bed body:	1 set



Optional

Castors:	optional
Color:	customizable



KL-TC003

Emergency Trolley



*Clapboard



KL-TC001

Emergency Trolley
Size:780*480*980mm



KL-TC002

Emergency Trolley
Size:680*525*1000mm



Technical Configuration

Main size:	750*480*920mm(small size)850*520*950(large size)	Guardrail:	1 pc
Main material:	ABS	Mute castors:	4 (2 brakes)
Working surface:	3 pcs	The storage space can be customized:	large/medium/small/etc.
6*4 swivel hook infusion stand:	1 set	Defibrillation plate,Basket, Dirt bucket, etc.:	optional

KL-TC004

Emergency Trolley
Size:750*475*930mm
850*520*960mm



KL-TC005

Anesthesia Trolley
Size:750*475*930mm
850*520*960mm



KL-TC006

Clinic Trolley
Size:750*475*920mm



KL-TC007

Medicine Trolley
Size:750*480*920mm
850*520*950mm



KL-TC020

Emergency Trolley
Size:750*475*930mm



KL-TC021

Anesthesia Trolley
Size:750*475*930mm



KL-TC022

Treatment Trolley
Size:750*475*930mm



KL-TC023

Emergency Trolley
Size:850*520*990mm



KL-TC008

Emergency Trolley
Size:625*475*930mm



KL-TC009

Anesthesia Trolley
Size:625*475*930mm



KL-TC010

Clinic Trolley
Size:625*475*930mm



KL-TC011

Medicine Trolley
Size:625*475*920mm



KL-TC024

Anesthesia Trolley
Size:670*450*880-960mm



KL-TC025

Clinical Trolley
Size:670*450*880-960mm



KL-TC026

Clinical Trolley
Size:670*450*825-1000mm



KL-TC027

Clinical Trolley
Size:670*450*880-960mm



KL-TC012

Case History Trolley
Size:450*410*950mm



KL-TC013

Case History Trolley
Size:750*410*950mm



KL-TC014

Case History Trolley
Size:450*492*1040mm



KL-TC015

Case clip
Size:230*330mm
Color can be customized



KL-TC028

Clinical Trolley
Size:775*635*910mm



KL-TC029

Clinical Trolley
Size:775*635*910mm



KL-TC030

Clinical Trolley
Size:780*475*1000mm



KL-TC031

Clinical Trolley
Size:775*630*1080mm



KL-TC016

ABS Trolley
Size:730*480*900mm



KL-TC017

ABS Trolley
Size:730*480*900mm



KL-TC018

ABS Trolley
Size:539*366*845mm
640*456*845mm
760*520*845mm



KL-TC019

ABS Trolley
Size:539*366*855mm



KL-TC032

Laundry Trolley
Size:800*600*840mm



KL-TC033

Laundry Trolley
Size:900*450*1050mm



KL-TC034

Laundry Trolley
Size:φ500X500X800mm



KL-TC035

Wast Collecting Trolley
Size:800*590*890mm



KL-TC036

Instrument Trolley
(Disassembly Structure)
Size:670*450*900mm



KL-TC037

Instrument Trolley
Size:670*450*900mm



KL-TC038

Instrument Trolley
Size L:1550*460*900mm



KL-TC039

Stainless Steel Trolley
Size L:1000*500*900mm
M: 880*430*860mm
S : 760*360*820mm



KL-TC052

Swab Stand
Size:900*500*1500mm



KL-TC053

Medicine Trolley
Size:670*450*1000mm



KL-TC054

Infusion Treatment Trolley
Size:670*450*1700mm



KL-TC055

Wet Linen Trolley
Size:815*610*750mm



KL-TC040

Dressing Trolley
Size:670*450*920-900mm



KL-TC041

Treatment Trolley
Size:460*460*860mm



KL-TC042

Treatment Trolley
Size:670*450*900mm



KL-TC043

Dressing Trolley
Size:670*450*900mm



KL-TC056

Shelf with Four Layers
Size:1200*600*1800mm



KL-TC057

Oxygen Bottle Trolley
Size:1435*360*1040mm



KL-TC058

Work Table
Size:1200*800*1600mm



KL-TC059

Work Table
Size:1200*800/300*1600mm



KL-TC044

Anesthesia Trolley
Size:600*400*850mm



KL-TC045

Insulated Food Cart
Size:1030*550*980mm



KL-TC046

Emergency Trolley
Size:700*400*850mm



KL-TC047

Emergency Trolley
Size:700*400*850mm




KL-TC060


Stainless Steel Washing Sink Size
1600(W)x630(D)x1350(H)mm

Technical Configuration & Optional


Faucet	2 pcs
Soap Dispenser	1 pc



Size:1000(W)x630(D)x1350(H)mm



Size:2000(W)x630(D)x1350(H)mm



Size:2400(W)x630(D)x590(H)mm



KL-TC061

Cabinet



KL-TC062

Cabinet



KL-TC048

Case History Trolley
(20 shelves)
Size:340*400*900mm



KL-TC049

Case History Trolley
(40 shelves)
Size:630*430*900mm



KL-TC050

Wash Basin
Size:Φ520*810-1300mm



KL-TC051

Mayo Trolley
Size:500*350*870-1300mm



KL-TC209

Ultrasonic cart



Please contact customer service for more customization

Technical Parameter

Top surface :	40.5*38cm
Basket :	30*17*15cm
Base size :	47*47cm
Wheel diameter :	7.5cm
Inner diameter of drawer :	30*24cm
Lower drawer support plate:	39.5x28cm
* Table-board can choose three size: 40.5*38 cm / 47*42 cm / 50*50 cm	



Technical Parameter & Optional

Cold rolled steel plate base:	1 set	Castors:	4 pcs
ABS dust cover medical silent pully:	4 pcs		
Gold bending basket:	2 pcs		
ABS flat with handle support plate:	1 set		
Support plate with drawer:	1 set		

KL-TC210

ECG cart



Technical Parameter

Top surface :	420*400mm
Base diameter:	500mm
Basket size:	320*200*120mm
Height:	830mm

Technical Parameter & Optional

Cold rolled steel plate base:	1 set
Medical silent pully:	4 pcs
ABS flat with handle support plate:	1 set

KL-TC211

Monitor cart

Technical Parameter

Table size:	32.8*33.5cm
Column diameter:	4.5cm
Wire basket size:	27.5*15.5*10cm
Knob can adjust the height, adjust the range:	145cm-185cm
Column:	Aviation aluminum, surface anode
Strap handle :	1 set
Pentagram base:	4 pcs
Universal wheel :	5



KL-TC201

Simple Laptop Cart

Size:535*45*750-1150mm



KL-TC202

All-in-One Computer Cart

Size:490*780*980-1030mm



KL-TC203

Monitor Cart

Size:510*1000mm



KL-TC204

Table Cart

Size:510*830mm



KL-TC205

Endoscopic System
Equipment Cart

Size:665*630*1650mm



KL-TC206

Information Nursing Cart

Size:520*620*1315mm



KL-TC207

Ward Rounds Cart

Size:530*520*730mm



KL-TC208

Patient Monitor Cart

Size:505*630*950mm





KL-MS001

Foldable Stretcher
Size:187x50x18cm



KL-MS003

Foldable Stretcher
Size:185x48x21cm



KL-MS002

Foldable Stretcher
Size:188x50x20cm



KL-MS004

Foldable Stretcher
Size:220*55*15cm



KL-MS005

Aluminum alloy plate stretcher
Size:185x50x3cm



KL-MS007

Plastic shovel stretcher
Size:162x44.5x7cm



KL-MS006

Scoop Stretcher
Size:210x44x6cm



KL-MS008

Hanging basket stretcher
Size:216x61x19cm



KL-MS009

Foldable Stretcher
Size:125x55x17cm



KL-MS011

Foldable Stretcher
Size:138x18x17cm
123x18x17cm



KL-MS010

Electric Stair Stretcher
Size:123x50x16cm



KL-MS012

Stair Stretcher
Size:215x54.5x13cm



KL-MS013

Exticaton Device

Packing Size:187x50x18cm



KL-MS015

Cervical Collar



KL-MS017

Spine Board



Size:1380*350*50mm

Size:1830*450*50mm

Size:2210*630*170mm

Size:1830*450*50mm

Size:1840*450*55mm

Size:1860*410*50mm

KL-MS014

Head Immobilizer



KL-MS016

Patient Transfer

Size:180x50x1cm



KL-MS018

Neil-Roberison Stretcher

Load Bearing:120kg



KL-MS019

Vacuumm Stretcher

Packing size:80x88x80cm



KL-MS020

Ambulance Stretcher



Technical Configuration & Optional

High Position (L*W*H):	190*55*92cm	Load Baring:	≤159kg
Low Position (L*W*H):	195*55*26cm	N.W:	39kg
Packing size:	195*62*28cm	G.W:	42kg
Maximum Angle:	75°	Color:	customizable

KL-MS022

Chair stretcher



Technical Configuration & Optional

High Position (L*W*H):	196*55*86cm	Castors:	4 pcs
Low Position (L*W*H):	196*55*25cn	Color:	customizable
Packing size:	198*62*26cm		
Maximum Angle:	85°		

KL-MS021

Automatic Loading Ambulance



Technical Configuration & Optional

Product size:	190*55*50cm	N.W:	29kg
Packing size:	195*62*34cm	G.W:	35kg
Maximum Angle:	75°	color:	customizable
Load Baring:	≤159kg		

KL-MS023

Automatic Loading Ambulance




Technical Configuration & Optional

High Position (L*W*H):	190*55*92cm	Load Baring:	≤159kg
Low Position (L*W*H):	195*55*26cm	N.W:	34.5kg
Adjusted size:	195*62*28cm	G.W:	40.5kg
Packing size:	195*62*42cm	Color:	customizable

KL-PT001

Patient Transportation Trolley



Technical Parameter

Length:	1960mm
Width:	600mm
Height:	750cm
Weight capacity:	150kg
Metal Castor wheels:	12.5cm
Stainless steel bed frame material:	1 set
Mattress with Learther cover high:	4cm



Technical Configuration & Optional

Guardrail :	1 set	Color of mattress:	customizable
IV pole with 2 hooks :	1 set		
Separation stretcher:	2 sets		
Castors:	4 pcs		

KL-PT002

Patient Transportation Trolley



Technical Parameter

Length:	2003mm
Width:	640mm
Height adjustment:	530-900mm
Back panel lifting:	0-85°
Safety load:	160kg
Head and tail size:	740*600*45mm
Mattress size:	1950*600*60mm



*Safety lock

Technical Configuration & Optional

Silent castors:	4 pcs	Castors:	4 pcs
Imported hydraulic cylinder:	1 set	Color:	customizable
Guardrail:	1 set		
PP bed panel:	3 sections		

KL-PT003

Patient Transportation Trolley



Technical Parameter

Length:	1930mm
Width:	640mm
Height adjustment:	590*870mm
Backrest up:	0-75°
Load capacity:	250 kg
Caster diameter:	6-inch
*Please refer to the actual object	



Technical Configuration & Optional

Adjustable IV rod:	1 pc	Silent castors:	4 pcs
Patient safety belt:	2 pcs	Crank system:	1 set
Oxygen cylinder rack:	1 set	Backplane air spring:	1 set
ABS guardrail:	1 set	Mattress:	1 pc (customizable)

KL-PT004

Patient Transportation Trolley



Technical Parameter

Length:	1930mm
Width:	660mm
Height adjustment:	510-850mm
Backpanel lifting:	0-70°(±10°)
*Please refer to the actual object	



Technical Configuration & Optional

Shaking hand system:	1 set	Holes for infusion stand:	2 pcs
Controllable gas spring:	1 set	Infusion stand:	1 set
Multifunctional PP guardrail:	1 set	Oxygen cylinder holder:	1 set
6-inch centralcontrol castors:	4 pcs	Central control brake pedal:	4 pcs
Guide castors:	1 set	Transfer mattress:	1 pc
Guide pedal:	2 pcs	Safety lock:	2 pcs

KL-PT005

Patient Transportation Trolley



Technical Parameter

Length:	2100mm
Width:	750mm
Height adjustment:	600-920mm
Backpanel lifting:	0°-80°
Leg lifting:	0°-30°
Reverse Trendelenburg position:	0-18°
Trendelenburg position:	0-18°
Mattress:	1 pc (customizable)
*Please refer to the actual object	



KL-PT006

Operation Connecting Trolley



Technical Parameter

Length:	3650mm
Width:	640mm
Height:	640-970mm
Backpanellifting:	0-70°
IV Pole:	1 Set
Mattress:	1 pc (customizable)
*Please refer to the actual object	



KL-LY001

Wheelchair

Technical Parameter

Lengh:	94cm
Height:	88cm
Folding width:	22cm
Unfolding width:	58cm
Seat width:	46cm
Seat height:	40cm
Wheel:	8"



KL-LY002

Wheelchair

Technical Parameter

Product size:	103*68*90cm
Folding width:	27cm
Seat Lengh:	42cm
Seat width:	45cm
Seat height:	50cm
Guard plate:	PP plastic guard plate
Handrail:	ABS hand pad



KL-LY003

Wheelchair

Technical Parameter

Product size:	103*68*85cm
Folding width:	27cm
Seat Lengh:	42cm
Seat width:	45cm
Seat height:	45cm
Load:	100KG
Brake:	Hand brake



KL-LY004

Electric Wheelchair



Technical Parameter

External (L*W*H):	1100*650*90mm
Sitting (H*W*D):	500*460*430mm
Speed:	6 km/h
Load capacity:	100kg
Motor:	2*250w
N.W.:	43KG
G.W.:	45KG



Technical Configuration & Optional

Ront wheel size:	10 inches	Color:	customizable
Rear wheel size:	16 inches		Battery optional
Driving distance:	15 km		
Battery:	lead-acid 12A		

KL-LY006

Wheelchair



Technical Parameter

Product size:	103*68*89cm
Folded size:	103*29*89cm
Cushion size:	42*46cm
Wheel:	7/24 inches
Wheelchair load capacity:	100kg
Wheelchair weight:	22KG
Adjustable backrest:	90-180



Technical Configuration & Optional

Wheelchair:	1 pc	Rear storage bag:	1 set
Commode:	1 set	More specifications can be consulted.	
Removable cushion:	1 set		
Leg rests:	2 pc		

KL-LY005

Electric Wheelchair



Technical Parameter

Product size:	850*440*790mm
Product weight:	17kg
Motor power:	250W*2
Battery capacity:	24V12A
Gradeability:	≤25°
Battery life:	about 12-15 kilometers
Wheel size:	10/16 inches



Technical Configuration & Optional

Wheel chair:	1 pc	Detachable leg cloth:	1 pc
360°controller:	1 set	Non-slip pedals:	2 pcs
Thicker tires:	2 pcs	Charger:	1 set
Universal wheels:	1 set		Battery optional

KL-LY007

Wheelchair



Technical Parameter

Product size:	103*49*92cm
Seat cushion (L*W):	37*37cm
Sitting height:	46cm
Load-bearing:	100kg
G.W.:	14kg
N.W.:	16kg



Technical Configuration & Optional

Side baffle:	2 pcs	Adjustable plastic dining table:	1 set
Adjustable handrail:	1 pc	Front wheel:	7 inches
Foot rest:	2 pcs	Rear-wheel:	22 inches
Seat belts:	1 set	Color:	customizable

KL-LY790

Bath Chair



Technical Parameter

Product height:	700-850mm
Seat height:	370-520mm
Seat depth:	380mm
Bearing capacity:	150KG
G.W.:	3kg
N.W.:	2.5kg

KL-FS792

Commode



Technical Parameter

Product size:	280*140*730-830mm
Seat height:	400-500mm
Seat depth:	400mm
Seat width:	450mm
Adjustable height:	100mm
Bearing capacity:	100KG
Net weight:	2.1 kg

KL-FS901

Armpit Crutch



Technical Configuration & Optional

Product size:	103-123cm/113-133cm/123-143cm	Bearing capacity:	150kg
Adapt to height:	130-153cm/154-174cm/175-195cm	Adjustment gear:	9 gears can be adjusted
Product material:	stainless steel	Quantity of crutches:	2 sticks
Tube wall thickness and width:	0.7mm	Detachable armrest:	2 sets
Product weight:	0.8kg	Anti-slip mat:	2 pcs

KL-FS902

Elbow Crutch



Technical Configuration & Optional

Packing size:	14cm*83cm*18.5cm	Load Weight :	150kg
Fold type:	17cm*45cm*15.5cm	Net weight(Retractable) :	0.45kg
Product material:	Aluminum elbow	Net weight(folding) :	0.48kg
Gear adjustment:	10 gears can be adjusted	Applicable height:	140-190cm
Tube wall thickness:	12MM	Color:	Silver/black paint

KL-FS903

Tri-fold Crutch



Technical Configuration & Optional

Product length:	123cm	Color:	customizable
Fold length:	30cm		
Product material:	stainless steel		
Tube wall thickness and width:	0.7mm		
Product weight:	0.3kg		

KL-FS904

Stainless Steel Crutch with Chair



Technical Configuration & Optional

Product size:	28*21*85cm	Open the bench height:	48cm
product weight:	About 0.5kg	Product color:	Blue/Brown/Black
Height of fold:	85cm		
Product bearing capacity:	130kg		
Product Material:	Stainless steel		

KL-FS905

Walking Aid



Technical Configuration & Optional

Product size:	50*43*75-87cm	Anti-slip mat:	2 pcs
Product material:	stainless steel	Castor:	2 pcs
Tube wall thickness:	1.2mm	The sitting board is removable	
Bearing capacity:	100KG		
Castor:	7 inch		

KL-FS907

Cart Walker



Technical Configuration & Optional

Product size:	68*61*80-90cm	Load-bearing:	120kg
Basket size:	40*20*15cm	N.W:	7.5kg
Size of seat cushion:	40*35cm	Wheel:	17.5cm
Sitting height:	50cm	Grip spacing:	47cm
Height of back:	35cm	Color:	customizable

KL-FS906

Walking Aid



Technical Configuration & Optional

Product size:	50*62*78-96cm
Folding Size:	55*10*78mm
Tube wall thickness :	1.2mm
Bearing capacity:	100KG
Adjustment gear:	6 gears can be adjusted

KL-FS908

Cart Walker



Technical Configuration & Optional

Product size:	73*40*91cm	Load-bearing:	120kg
Basket size:	40*40*30cm	N.W:	7.5kg
Size of seat cushion:	41.5*41.5cm	Wheels:	17.5/30cm
Sitting height:	48cm	Grip spacing:	47cm
Height of back:	42cm	Color:	customizable

KL-CE001

Certrifuge



Technical Parameters & Configuration

Maximum speed:	4000 RPM
Maximum centrifugal force:	2325XG
Timing:	0-30min
Rotor capacity:	20mlX12 (standard configuration)
Power supply:	110V 60HZ 220V 50Hz 110V
Machine size(mm):	280*310*265
Packing size (mm):	340*350*330
Gross/net weight:	10/9KG
Body:	1
Non-slip mats:	4
Aluminum alloy centrifugal rotor:	12 holes

KL-DT001

Electrothermal Constant-Temperature
Drying Oven



Technical Parameters & Configuration

Product size:	65*46*46CM
Inner container size:	35*35*35CM
Temperature range:	indoor temperature ±10℃ -250℃
Instrument accuracy:	1℃
Power supply voltage:	220V/50HZ
Heating power:	800W
Timing time:	0-999min
Liner material:	Galvanize
Standard partition:	Galvanized separator*2

KL-ST001

Portable Pressure Steam Sterilizer



Technical Parameters & Configuration

Volume:	45L
Power:	2.5KW
Voltage:	110v/220v 60Hz/50Hz
Rated working pressure:	0.25Mpa
Rated operating temperature:	139℃
Sterilization chamber volume:	316x615
Dimensions:	657*870*575
Packing size:	980*750*630
Gross weight/net weight:	142/105KG
Open top water tank:	1
Built-in printer:	1
Built-in steam generator:	1

KL-ST002

Vertical Pressure Steam Sterilizer



Technical Parameters & Configuration

Chamber volume:	35L ϕ318 x450mm
Working pressure:	0.22Mpa
Working temperature:	134℃
Adjustment of temperature:	40℃-134℃
Time for sterilization:	0~60min
Time for drying:	0~60min
Heat average:	≤+2℃
Power:	9kw/380V 50Hz
Dimension:	1400x600x1300mm
Transport dimension:	1550x750 x1850mm
G.W/N.W:	320/240kg
More specification consultation	

KL-ST003

Portable Pressure Steam Sterilizer



Technical Parameters & Configuration

Sterilizing volume:	18L ϕ280 x260mm
Working pressure:	0.14-0.16Mpa
Working temperature:	126℃
Adjustment of temperature:	105-126℃
Time for sterilization:	0~99min
Time for drying:	0.165Mpa
Heat average:	≤1℃
Power:	AC220V,50Hz/2KW
Transport dimension:	420*420*630
G.W/N.W:	14/12kg
More specification consultation	

KL-ST004

Horizontal Pressure Steam Sterilizer



Technical Parameters & Configuration

Chamber volume:	150L ϕ440×1000mm
Working pressure:	0.22Mpa
Working temperature:	134℃
Adjustment of temperature:	40℃-134℃
Time for sterilization:	0~60min
Time for drying:	0~60min
Heat average:	≤+2℃
Power:	9kw/380V 50Hz
Dimension:	1400x600x1300mm
Transport dimension:	1550x750 x1850mm
G.W/N.W:	320/240kg
More specification consultation	



WICKLAND MEDIKA IMANTARA

MEDICAL EQUIPMENT , AESTHETICS AND LABORATORY SUPPLIER



LOCATION

PALEM SEMEMI RAYA NO. 12-12A
KEL. SEMEMI, KEC. BENOWO, SURABAYA
JAWA TIMUR, INDONESIA



CONTACT US

031-33300312
031-33100412



EMAIL

WICKLANDMEDIKA@GMAIL.COM