

**WICKLAND MEDIKA IMANTARA**

MEDICAL EQUIPMENT , AESTHETICS AND LABORATORY SUPPLIER

CH



**CH-13 Series**  
FLEXIBLE VIDEO  
CHOLEDOCHOSCOPE

## CMOS+LED Technology

Fiber-free design makes it more durable.



## High-quality Image

Excellent image gives more details.



## Multi-function Buttons

Improve operation flexibility and accuracy.

- ① Zoom in/out
- ② Capture/record
- ③ Freeze/release



## Portable

Equiped with image processor supports external monitor.



## Specifications

Model	Insertion Tube	Distal End	Working Channel	Working Length	Field of View	Depth of Field	Bending (Up/Down)
CH-1331	3.1mm	3.2mm	1.2mm	450mm	120°	3-50mm	130°/130°
CH-1349	4.8mm	5.0mm	2.0mm	380mm	120°	3-50mm	175°/130°

## Image Processor



10.1" Touch screen



1280\*800 HD Resolution



HDMI+ AV output



Rechargeable battery supporting 4 working-hour



Detachable 64G SD card



### LOCATION

PALEM SEMEMI RAYA NO. 12-12A  
KEL. SEMEMI, KEC. BENOWO, SURABAYA  
JAWA TIMUR, INDONESIA



### CONTACT US

031-33300312  
031-33100412



### EMAIL

WICKLANDMEDIKA@GMAIL.COM