



WICKLAND MEDIKA IMANTARA

MEDICAL EQUIPMENT , AESTHETICS AND LABORATORY SUPPLIER



CHD-1 series

Single-Use Flexible
Cholangioscope

Dual light sources, dedicated irrigation,
large 2.0mm working channel.



Ergonomic and lightweight,
with deflection lock reduces fatigue.



Provides four-way steering,
provides traditional sens of operation.



Attaches to duodenoscope for single operator use.



Exclusive channels for different operations



Instrument,
Irrigation flow control



Suction controlled



Suction



Water jet /
Guide wire insertion



Reduces cost and risk associated with
reprocessing and cross-contamination.

Product Parameter

Model	Insertion Tube	Working Channel	Water jet Channel	Suction Channel	Working Length	Bending (U/D/L/R)	Field of View	Depth of Field
CHD-1	3.7mm	2.0mm	0.6mm	0.6mm	2120mm	120°	120°	3-50mm

Image Processor



HV-3101
Image Processor



10.1 inches
Touch screen



1280*800
High Resolution



HDMI+ AV
output



Detachable rechargeable battery
supporting 4 hours working



Detachable
SD card



Portable and
lightweight design



LOCATION
PALEM SEMEMI RAYA NO. 12-12A
KEL. SEMEMI, KEC. BENOWO, SURABAYA
JAWA TIMUR. INDONESIA



CONTACT US
031-33300312
031-33100412



EMAIL
WICKLANDMEDIKA@GMAIL.COM