

Ear, Nose and Throat Examination and Treatment Table

Technical Parameters and Characteristics

The YW-1001 Ear, Nose and Throat Examination and Treatment Table is manufactured using advanced technology

- The countertop is made of tempered glass material, and the body is made of cold rolled steel plate with surface spray molding technology
- Corrosion resistant and rust resistant, with long service life
- Adopting oil-free compressor and vacuum pump, low noise and long service life
- The maximum pressure of the positive pressure pump can reach 400Kpa, with automatic shutdown function
- Automatically work when the pressure drops below 400Kpa
- The maximum negative pressure of the vacuum pump can reach -0.09Mpa, with strong power
- This machine is equipped with 4 spray guns, 1 suction gun, and 1 blowing gun
- Convenient for doctors to use treatment flexibly
- Spray gun and blowing gun pressure: 100Kpa-150Kpa adjustable
- Suction gun negative pressure: 0-0.07Mpa adjustable
- Light source illumination: >1000Lux

YW-1001 configuration

- 1 main unit (length 1.7 meters with writing desk)
- 4 spray guns (2 straight and 2 curved)
- 1 suction gun and 1 blowing gun
- 1 set of mist removal device
- 1 stainless steel tray
- 2 stainless steel disinfection cylinders
- 4 glass reagent bottles (2 white, 2 brown)
- 12500ml dirt bottle

Optional: camera system, graphic workstation, color printer, otoscope, nose mirror, laryngoscope



YW-1002 parameters are the same as YW-1001

YW-1002 configuration

- 1 main unit (length 730mm)
- 4 spray guns (2 straight and 2 curved)
- 1 suction gun, 1 blowing gun
- 1 set of mist removal device
- 2 stainless steel disinfection cylinders
- 4 glass reagent bottles (2 white and 2 brown)
- 1 2500ml dirt bottle

Optional:

- Camera system
- Image and text workstation
- Color printer
- Far glasses, nose glasses, and laryngoscope

