

EMK-A1

FULL HD Endoscope System



SPEIFICATION

Camera head	• Full HD imaging (1920 x 1080) • Full HD optical coupler C-mount standard
Frequency	• 50/60 Hz
Device protection	• Emi Filter, 2x Fuse, 3A / 250V
Video system	• Image resolution: FHD 1920 x 1080 • Scanning system: Progressive Scanning • White balance automatically
Screen size	• 5 inch
Recording (optional)	• Image format: JPG 1920 x 1080 • Video format: Full HD 1920 x 1080 • MPEG-4 (H.264 Codec); Image rate 60fps • Connection: USB
Connections	Video out terminal • 1x HDMI
Protection class	• I
Safety class	• BF
Standard	• IEC 60601-1:2014 • IEC 60601-2-18
Features	• White balance • Freeze function • Zoom in / Zoom out • Touch Buttons / Touch Screen • USB storage • Photo capture / Video recording
Language	• Chinese/English/Japanese/Korean/Spanish/German/Russian
Dimensions	• 356(L) x 366(W) x 135(H)
Weight	• 7,3 Kg