


HYSTEROSCOPY



PT. Wickland Medika Imantara
Medical, Aesthetics & Laboratory Supplier
Jl. Palem Sememi Raya No 12 - 12A Kel.
Sememi,
Kec. Benowo Surabaya 60198 | 031-33100412

Hysteroscope (Z-Shaped) SWD-GQJ-II




Hysteroscope SWD-GQJ-VI

Y-shaped Cold knife

GQJ-VI.30.07

fit $\phi 3 \times 360\text{mm}$ instruments



Hysteroscope


(Mirena Suture Surgical Instruments Set)



Optional Instruments ($\phi 4 \times 360\text{mm}$)

	hook scissors	208.03.04		monopolar electric coagulator	235.01.04
	curved scissors	209.04.04		monopolar electric hook	234.01.04
	straight scissors (single action)	210.03.04		monopolar electric hook	234.02.04
	straight scissors (double action)	210.04.04		bipolar electric coagulator	235.03.04
	dissecting forceps	213.01.04		bipolar electric hook	234.03.04
	Maxi forceps	220.09.04		bipolar electric spud	236.03.04
	biopsy forceps	220.14.04		high frequency cable	270.01
	claw forceps 2+3 teeth	221.12.04		high frequency cable	270.02

Hysteroscope (integrated)



Hysteroscope(integrated) SWD-GQJ-III




Item No.	Name	Visual direction	Working length	Diameter	Instrument channel
GQJ30.02	Hysteroscope	30°	205mm	6mm	5Fr, 1.6x400mm
GQJ30.03				6.4mm	7Fr, 2.0x400mm

semirigid instruments SWD-GQJ (360° rotatable, detachable into handle and insert)

Disassembled

Disassembled in 4 steps ① unscrew the rotated wheel ② open the handle ③ ④ pull out the insert up-turned



φ1.6x400mm(5Fr blue)	φ2.0x400mm(7Fr green)	note
foreign body forceps G2301.05	foreign body forceps G2301.07	
scissors G2302.05	scissors G2302.07	single action, blunt
biopsy forceps G2303.05	biopsy forceps G2303.07	
grasper G2305.05	serrated forceps G2305.07	
scissors G2306.05	scissors G2306.07	single action, pointed
scissors G2304.05	scissors G2304.07	double action (blunt)
foreign body forceps 2307.05	/	
electric coagulator III B235.01.16		monopolar
electric coagulator IV B235.04.16		bipolar
non-detachable instrument		optional

non-detachable instrument



Hysteroscope



Hysteroscope GQJ30.29
2.9x302mm 30°



Hysteroscope GQJ30.40
4x302mm 30°



Manipulator GCZQ01.01
5.4x200mm



Manipulator GCZQ02.01
6.4x200mm



Sheath&obturator GJQ01.02
6.2x185mm



Sheath&obturator GJQ02.02
7.4x185mm

Surgical instruments

The size of instruments is 5Fr and can be detached into handle and insert



Hysteroscope (SWD-GQJ)



GQJ22.64 22°



GQJ00.64 0° / GQJ30.64 30°

Disassembled

Disassembled in 4 steps ①unscrew the rotated wheel ②open the handle ③④pull out the insert up-turned



flexible biopsy
forceps A2101



flexible foreign body
forceps A2102



flexible grasper
A2105



flexible flexible serrated
forceps A2103



flexible scissors
A2104

cup-shaped uterine manipulator SWD-F

Type-A F700.01



F703.01 measuring bar



F703.01.01

F703.01.02



F703.01.03

F703.01.04



F706.01 cone head



F705.01 spiral myoma drill



F705.02 spiral myoma drill



F704.01 booster rod



F702.01 manipulator rod



F707.01 Cup(L)



F707.02 Cup(M)



F707.03 Cup(S)



F701.01 holder

cup-shaped uterine manipulator SWD-F

Type-B F 700.02



F703.02 measuring bar



F706.02 cone head



F705.01 spiral myoma drill



F705.02 spiral myoma drill



F704.01 booster rod



F702.01 manipulator rod



F707.01 Cup(L)



F707.02 Cup(M)



F707.03 Cup(S)



F701.01 holder

uterine cannula with irrigation SWD-E

E425.02 uterine cannula with irrigation length: 400mm



liquid channel

cone head

cone pin

slider

valve