

Endoscopic surgical shaver system



PT. Wickland Medika Imantara

Medical, Aesthetics & Laboratory Supplier

Jl. Palem Sememi Raya No 12 - 12A Kel. Sememi,
Kec. Benowo Surabaya 60198 | 031-33100412

SPECIFICATION

PRODUCT INTRODUCTION

Application: Used in endoscopic surgery to achieve the surgical function of suction and excising soft tissue.

Indications: Could be used with laser and resectoscopes.

FUNCTION CHARACTERISTIC

RPM: 1000-8000 r/min adjustable

Oscillating rotation mode available

Suction power adjustable

7 inches touch screen

Sterilization: Autoclave or low-temperature plasma

Waterproof handle and footswitch

Bottle with overflow valve

TECHNICAL PARAMETERS

Voltage: 110V - 220V 50/60Hz

RPM: 1000 - 8000 r/min, 1000 per step

Aspiration: -20 ~ -60 kPa, 1 kPa per step

Maximum torque: not less than 60mN.M.

Blade: working length 365mm, ϕ 4.5mm

ACCESSORIES

Blade (Customization available)



Handle



Foot switch