



ENDOUROLOGY UPPER LOWER




PT. Wickland Medika Imantara

Medical, Aesthetics & Laboratory Supplier

Jl. Palem Sememi Raya No 12 - 12A Kel. Sememi,
Kec. Benowo Surabaya 60198 | 031-33100412

DAFTAR ISI	
Daftar Nama	Halaman
Cystoscopy	1
Resectoscope	2-3
Uretero-Renoscropy	4-6
Nephroscope	7
Mini-Neproscope	7
Monopolar Resectoscope	8-9

CYSTOSCOPE



Cystoscope

PGJ00.01 $\Phi 4 \times 302\text{mm}$ 0°
PGJ30.01 $\Phi 4 \times 302\text{mm}$ 30°
PGJ70.01 $\Phi 4 \times 302\text{mm}$ 70°

Disassembled

Disassembled in four steps ①unscrew the rotated wheel ②open the handle ③④pull out the insert up-turned



RESECTOSCOPE

Monopolar Resectoscope

SWD-DQJ-I



- outer sheath DQJ01.02.29



- outer sheath DQJ01.02.40



- endoscope DQJ301.29 $\phi 4 \times 302 \text{mm}$ 30°



- working element
DQJ01.01.29



- inner sheath
DQJ01.03.29



- obturator
DQJ01.04.29



Monopolar



4030
straight loop



4023
roller-shaped



4028
ball-shaped



4026
needle-shaped



4021
cutting loop



4022
spud-shaped



4024
spike-shaped

- electric vacuum
DQJ01.08.01



- adaptor
DQJ01.07.01



- high frequency cord
DQJ01.06.29



URETERO-RENOSCOPY

Ureteroscope SWD-SNJ-I



product name	item NO.	specification	working length	instrument channel	Fr	note
ureteroscope	SNJ-I-003430	3mm	430mm	3Fr(1.0mm)	4.5/6.5Fr	seminal cystoscope
ureteroscope	SNJ-I-003315	3mm	315mm	3Fr(1.0mm)	4.5/6.5Fr	seminal cystoscope
ureteroscope	SNJ-I-003250	3mm	250mm	3Fr(1.0mm)	4.5/6.5Fr	seminal cystoscope

Ureteroscope SWD-SNJ-I



product name	item NO.	specification	working length	instrument channel	Fr	note
ureteroscope	SNJ-I-004430	4mm	430mm	4Fr(1.3mm)	6/7.5Fr	Pediatric
ureteroscope	SNJ-I-004315	4mm	315mm	4Fr(1.3mm)	6/7.5Fr	Pediatric



single channel

Ureteroscope SWD-SNJ-I



product name	item NO.	specification	working length	instrument channel	Fr	note
ureteroscope	SNJ-I-005430	4mm	430mm	5Fr(1.6mm)	8/9.8Fr	adult
ureteroscope	SNJ-I-005315	4mm	315mm	5Fr(1.6mm)	8/9.8Fr	adult




single channel

Ureteroscope SWD-SNJ-I




Double channel ureteroscope,
Independent water in and out channels to reduce the risk of reverse infection.



product name	item NO.	specification	working length	instrument channel	Fr	note
ureteroscope	SNJ-I-004430S	4mm	430mm	4Fr(1.3mm)	6/7.5Fr	double channel



double channel



product name	item NO.	specification	working length	instrument channel	Fr	note
ureteroscope	SNJ-I-005430S	4mm	430mm	5Fr(1.6mm)	8/9.8Fr	double channel




double channel

Disassembled

Disassembled in 4 steps

①unscrew the rotated wheel②open the handle③④pull out the insert up-turned



non-detachable



product name	item NO.	specification	working length	wheel color	Fr	note
foreign body forceps	A2201.10	φ1.0mm	580mm	non-detachable	3Fr	4.5/6.5Fr (seminal cystoscope)
	A2201.13	φ1.3mm	580mm	purple	4Fr	6/7.5Fr(pediatric)
	A2201.16	φ1.6mm	580mm	golden	5Fr	8/9.8Fr(adult)
serrated forceps	A220.10	φ1.0mm	580mm	non-detachable	3Fr	4.5/6.5Fr (seminal cystoscope)
	A2202.13	φ1.3mm	580mm	purple	4Fr	6/7.5Fr(pediatric)
	A2202.16	φ1.6mm	580mm	golden	5Fr	8/9.8Fr(adult)
biopsy forceps	A2203.10	φ1.0mm	580mm	non-detachable	3Fr	4.5/6.5Fr (seminal cystoscope)
	A2203.13	φ1.3mm	580mm	purple	4Fr	6/7.5Fr(pediatric)
scissors	A2203.16	φ1.6mm	580mm	golden	5Fr	8/9.8Fr(adult)
	A2204.13	φ1.3mm	580mm	purple	4Fr	6/7.5Fr(pediatric)
	A2204.16	φ1.6mm	580mm	golden	5Fr	8/9.8Fr(adult)

NEPHROSCOPE

nephroscope
SWD-SNJ-I



product name	item NO.	specification	working length	instrument channel	Fr	note
nephroscope	SNJ-I-006250	4mm	250mm	6Fr(2.0mm)	8.5/12Fr	percutaneous nephroscope

MINI-NEPHROSCOPE

mini-nephroscope



Product name	Working length	Instrument	Fr	Note
mini-nephroscope	224mm	3.5mm	18Fr, 0°	percutaneous nephroscope

Monopolar Resectoscope



PT. Wickland Medika Imantara
Medical, Aesthetics & Laboratory Supplier

Jl. Palem Sememi Raya No 12 - 12A Kel. Sememi,
Kec. Benowo Surabaya 60198 | 031-33100412

Monopolar Resectoscope

SWD-DQJ-I



- outer sheath DQJ01.02.29



- outer sheath DQJ01.02.40



- endoscope DQJ301.29 $\phi 4 \times 302 \text{mm}$ 30°



- working element
DQJ01.01.29



- inner sheath
DQJ01.03.29



- obturator
DQJ01.04.29



Monopolar



4030
straight loop



4023
roller-shaped



4028
ball-shaped



4026
needle-shaped



4021
cutting loop



4022
spud-shaped



4024
spike-shaped

- electric vacuum
DQJ01.08.01



- adaptor
DQJ01.07.01



- high frequency cord
DQJ01.06.29

