

UMK-URS1

Uroflowmetry



Product Features

- Touch screen
- Support compatible printer
- Auto calibration
- Clean user interface
- Realtime display

Technical Specification

Voltage range (current consumption)	• 210-240 VAC
Frequency	• 50 Hz
Rated Current	• 0,6 A
Body Type	• ABS Plastic Body
Dimension	• 250mm x 250mm x 110mm (L x B x H)
Power Consumption	• 30 W
Temperature	• 5-40 °C
Sensor	• 0,5 mL weight system
	• 7" Touchscreen LCD Display for operating system
Accessories	• 1 Unit Uroflowmetry System • 1 Unit 7" Touch LCD Display • 1 Unit Keyboard & Mouse • 1 Unit Sensor & Measuring Cup • 1 Pc Pen Drive SanDisk • 1 Unit Printer

