



Wickland



CHD-1 series

Single-Use Flexible
Choledochoscope

Dual light sources, dedicated irrigation,
large 2.0mm working channel.



Provides four-way steering,
provides traditional sens of operation.



Ergonomic and lightweight,
with deflection lock reduces fatigue.



Attaches to duodenoscope for single operator use.



Exclusive channels for different operations



Instrument,
Irrigation flow control



Suction controlled



Suction



Water jet /
Guide wire insertion



Reduces cost and risk associated with
reprocessing and cross-contamination.

Product Parameter

Model	Insertion Tube	Working Channel	Water jet Channel	Suction Channel	Working Length	Bending (U/D/L/R)	Field of View	Depth of Field
CHD-1	3.7mm	2.0mm	0.6mm	0.6mm	2120mm	120°	120°	3-50mm

Image Processor



◀ HV-3101
Image Processor



10.1 inches
Touch screen



1280*800
High Resolution



HDMI+ AV
output



Detachable rechargeable battery
supporting 4 hours working



Detachable
SD card



Portable and
lightweight design