

Medical Endoscopy System Camera (Affordable Model With Recorder)

TU-Y5D100B (FHD)



Camera head

Image Sensor	SONY CMOS
Shortcut Key Function	Freeze/White balance/Digital zoom /Custom
FREEZE	Support
AWB(White Balance)	Support
D-ZOOM	0-8x
Working Environment Temperature	0°C~40°C

Endoscope Camera System

Input Interface	1080P
Output Interface	HDMI/ SDI/DVI/VGA/CVBS
U Disk	U disk USB2.0/USB3.0 (up to 32GB) to record and take photos
Power Supply	AC110V~240V
Resolution	4K(3840×2160 pixels) 1080P(1920×1080 pixels)
Brightness	1-30 steps
Support Lens Mode	PTED/Fiber/Arthroscope/Cystoscope/Hysteroscope/Laparoscope/ENT/Custom
Support Language	Chinese/English/Japanese/Korean/Etc
SNR(Signal Noise Ratio)	<45dB
Working Environment Temperature	0°C~40°C
Size	L320*W235*H90mm
Net Weight	3.5KG