



USB Output

*Compatible to hardware with USB video input. Everywhere E.N.T.



* EN-E19H Only

Portable

*Compatible with portable image processor.



* EN-E19H incompatible

Multi-function buttons


- ① Zoom in/out (*White balance)
- ② Capture/record (*Brightness)
- ③ Freeze/release



* EN-E19H Only

CMOS+LED Technology

CMOS image sensor delivers a high quality image.
No fiber: durable and easy maintenance.



Specifications

Model	Insertion Tube	Distal End	Working Channel	Working Length	Field of View	Depth of Field	Bending (Up/Down)
EN-E19H	3.4mm	3.6mm	\	300mm	100°	3-100mm	140°/140°
EN-1029	2.2mm	2.5mm	\	300mm	90°	3-50mm	130°/130°
EN-1039	3.1mm	3.4mm	\	300mm	90°	3-50mm	130°/130°
EN-1249	4.8mm	5.2mm	2.0mm	365mm	90°	3-50mm	130°/130°

Image Processor



10.1"
Touch screen



1280*800
HD resolution



HDMI+AV
output



Rechargeable battery
supporting 4 working hours



Portable
design



Detachable 64GB
SD card



PT. Wickland Medika Imantara
Medical, Aesthetics & Laboratory Supplier
Jl. Palem Sememi Raya No 12 - 12A Kel. Sememi,
Kec. Benowo Surabaya 60198 | 031-33100412