




AF

 **AF Series**
MOBILE FLEXIBLE VIDEO
BRONCHOSCOPE

CE 0197

Angle lock

Lock or release the angle of distal end for accurate diagnosis and treatment



Suction valve retainer

Unique design prevents the suction valve from falling off



Rotation

Improve operation flexibility and accuracy



Multi-function buttons

- ① Zoom in/out
- ② Capture/record
- ③ Freeze/release



Specifications

Model	Insertion Tube	Distal End	Working Channel	Working Length	Field of View	Depth of Field	Bending (Up/Down)	Application
AF-2031	2.8mm	3.0mm	1.2mm	610mm	120°	3-200mm	210°/130°	Anesthesia Pediatric
AF-2042	3.9mm	4.1mm	2.0mm	610mm	120°	3-200mm	210°/130°	Pediatric
AF-2049	4.8mm	5.2mm	2.0mm	610mm	120°	3-200mm	210°/130°	Pneumology
AF-2052	5.5mm	5.4mm	2.8mm	610mm	120°	3-200mm	180°/130°	ICU Anesthesia
AF-2055	5.8mm	6.0mm	2.8mm	610mm	120°	3-200mm	180°/130°	Pneumology
AF-2058	5.8mm	6.2mm	3.0mm	610mm	120°	3-200mm	180°/130°	Pneumology

Image Processor



◀ HV-50

- 720*720
- Detachable &
- Rechargeable
- lithium battery



▼ HV-3101

- 1280*800
- Touch screen
- HDMI+AV output
- Rechargeable battery supporting 4 working-hour



Touch screen



120° Up/Down



64GB TF card

