

PORTABLE MEDICAL HOLMIUM LASER (HO:YAG)



PT. Wickland Medika Imantara
Medical, Aesthetics & Laboratory Supplier
Jl. Palem Sememi Raya No 12 - 12A Kel. Sememi,
Kec. Benowo Surabaya 60198 | 031-33100412

URETEROSCOPIC LITHOTRIPSY

- 200um ultra-soft fibre optic with 360 degree rotation of the soft lens.

CLINICAL USE

- Nephrolithotomy
- Ureterolithotomy
- Cystolithectomy

TECHNICAL PARAMETERS

- Maximum single pulse output energy: 0.5-3.0 J
- Average output at terminal: 30W
- Power supply: 220V/50Hz
- Maximum operating frequency: 1-15 Hz
- Pulse width: 300us/600us

STANDARD FITTING



CONTROL PANEL

