



WICKLAND SINGLE USE URETEROSCOPE

IMAGE PROCESSOR

1080P HIGH RESOLUTION CAMERA



- Support USB drive and hard drive.
- Support FAT 32 and NTFS Format
- Support Recourding Function

PORTABLE PROCESSOR

EASY TO BE CARRIED WITH YOU




- 12.1 inches high resolution screen.
- Only weight 2.2kg.
- Support built-in battery solution.



CAMERA

160,000 PIXELS HIGH RESOLUTION



Model	Channel O.D.	Working Channel O.D.	Working Length
MD-A-UE30	3.0mm (9.0Fr)	1.2mm (3.6FR)	670mm
MD-A-UE25	2.5mm (7.5Fr)	1.2mm (3.6FR)	670mm



PT. Wickland Medika Imantara

Medical, Aesthetics & Laboratory Supplier

Jl. Palem Sememi Raya No 12 - 12A Kel. Sememi,
Kec. Benowo Surabaya 60198 | 031-33100412