



KL-LED.D78/D78.WMI

WICKLAND SHADOWLESS OPERATING LAMP

KL-LED.D78/D78.WMI

KL-LED.D78/D78.WMI series LED shadowless operating lamp is suitable for various operation needs. It is a modern ideal surgical lighting equipment operating room and purify the operating room.



INTELLIGENT

Intelligent LCD touch control



ENERGY-SAVING

imported LED cold light source



DURABLE

Professional quality certification

SPECIFICATIONS

- (1) Lamp head diameter: (overall) $\geq 72(\text{cm})$
- (2) Illuminance at 1 meter of lamp head $160000 \geq E_c \geq 80000/160000 \geq E_c \geq 80000$
- (3) Spot diameter D10: $180 \pm 40(\text{mm})/180 \pm 40(\text{mm})$
- (4) Spot diameter D50: $110 \pm 35(\text{mm})/110 \pm 35(\text{mm})$
- (5) Illumination depth: $1500(\text{mm})$
- (6) The color development index Ra was: ≥ 96
- (7) Red indicates R9: $\geq 95(\text{Ra})$
- (8) Optional color temperature: $6700 \geq T_c \geq 3000$ (multi-level adjustable)
- (9) Lamp life 50,000(hours)
- (10) Total irradiance: $< 2(\text{mW}/\text{m}^2\text{lx})$
- (11) Power supply voltage: $220\text{V} \pm 20\%$, 50Hz
- (12) Input power: $300 \pm 10\%$ (VA)
- (13) Installation height: $270-320(\text{cm})$

BASIC CONFIG

- (1) Lamp body: 2 pieces
- (2) Fixed plate: 1 piece
- (3) Rotating body: 1 set
- (4) Balance arm: 2 pieces
- (5) Fixed disc cover: 1 set
- (6) Switching power supply: 2



■ Nano coating of heat sink has created an excellent heat dissipation effect.

■ Intelligent touch screen, with adjustable color temperature


■ Titanium alloy balance arm, no welding, never rust

■ Imported spring arm, lift up and down, light and flexible

■ Minimally invasive lighting and R9 function, endoscopy function

MAIN FEATURE

KL-LED.D78/D78.WMI



DIVERSIFIED CUSTOMIZED SERVICES

The configuration of shadowless lamps can be customized and selected according to user needs, such as:

- Double-headed or single-headed shadowless lamp
- External boom camera
- Boom monitor

MULTI-SCENARIO APPLICATION



GENERAL HOSPITAL



DENTAL HOSPITAL



PET HOSPITAL



COSMETIC HOSPITAL



MEDICAL TEACHING