



PT. Wickland Medika Imantara
Medical, Aesthetics & Laboratory Supplier

Jl. Palem Sememi Raya No 12 - 12A Kel. Sememi,
Kec. Benowo Surabaya 60198 | 031-33100412

VLS-H11 Series Video Laryngoscope



Video Laryngoscope

Assist fast oral intubation with endotracheal tube



Compatible with HV-3101
10.1\" monitor for various
application scenarios.

HV-3101
· 1280*800
· Touch screen
· HDMI+AV output
· Rechargeable battery
supporting 4 working-hour

► VLS-H11 Series Specifications


Model	Patient Type
VLS -H11S	Children
VLS -H11M	Adult
VLS -H11L	Large adult
VLS -H11SML	Children, adult and large adult

► Product Features



Multi-function

3 function buttons make it easier to operate without contamination risk of the monitor



1:1 full screen

Restore the real picture



High compatibility

One image processor and 3 blade sizes make a seamless switching between adult and children.



High security

Aviation socket design secures solid connection.



Ergonomic design

Comfortable, portable and durable



Simple to operate

Low learning cost



Anti-fog design

Power on to use, no preheating.

► Application: Rapid oral intubation, oral examination and treatment.

► Applicable departments: Anesthesiology, emergency department, ICU and ambulance etc.

Disposable blade

Prevent cross infection

Patent No.: ZL 2022 2 0708531.6

1

Anti-fog coating

2

Easy installation

3

Friendly design

4

High quality

ALC-H1BL

ALC-H1ML

ALC-H1SL

