

Single Use Flexible Bronchoscope

Atraumatic Access • High Clarity • Single-Use Convenience



220° Up/Down



Rotatable



600ml/min^①
Powerful Suction



Soft Uniform Light



Encrypted Connector



PT. Wickland Medika Imantara

Medical, Aesthetics & Laboratory Supplier

Jl. Palem Sememi Raya No 12 - 12A Kel. Sememi,
Kec. Benowo Surabaya 60198 | 031-33100412



Compatible Accessories



1.8mm Biopsy Forcep



1.8mm Grasping Forcep



1.8mm Polypectomy Snare



1.8mm Cytology Brush

*More Respiration products are in developing


SINGLE USE CHOLANGIOPANCREATOSCOPE PARAMETERS

Model	MD-A-BO3212	MD-A-BO4828
O.D. of Tip	3.2mm (9.6 Fr)	4.8mm (14.4 Fr)
Bending Angle	U220° /D220°	U220° /D220°
Articulation Section O.D.	3.3mm (9.9 Fr)	5.0 mm (15 Fr)
Length	600mm	600mm
Working Channel	1.2mm (3.6 Fr)	2.8mm (8.4 Fr)
Depth of Field	3mm - 50mm	3mm - 50mm
Field of View	120°	120°
Light Source	LED	LED
Pixels Output	160,000	160,000
Single Use	Yes	Yes


IMAGE PROCESSING SYSTEM

One-Click Calibration • Clear Visualization • User-Friendly


Auto-Exposure Compensation for Dark Areas




Automatic Bad Pixel Correction



Intelligent Vertical Stripe Elimination




Intelligent Image Quality Enhancement



PARAMETERS



MODEL		MD-A-IP3010
Size	300mm×350mm×100mm	
Weight	5.1kg	
Signal Output	CVBS	Support
	HDMI 1.4	Support
	3G-SDI	Support
	DVI	Support
Resolution Output	1980X1080	
Frames/Sec Output	50-60HZ	
Signal Input	Type-C (Encrypted)	
Storage	USB 2.0 Storage Input	
Languages Support	English, Chinese etc.	



Built-in Intelligent Chip Can Automatically Achieve Image Quality Enhancement.



Encrypted Interface: Compatible with MEDNOVA's Various Electronic Endoscopes for Multiple De-



One-Touch White Balance to Prevent the Occurrence of Color Deviation.



Multiple Optional Output Interfaces, Including 3G-SDI for Ultra-Long-Range Signal Transmission

MODEL

		MD-A-PP285247
Size	285mmX50mmX247mm	
Screen Size	12.1 inches	
Weight	2.5kg	
Signal Output	HDMI 1.4	
Resolution Output	1980X1080	
Frames/Sec Output	50-60HZ	
Signal Input	Type-C (Encrypted)	
Storage	USB 2.0 Storage Input	
Languages Support	English, Chinese etc.	
Electronic Amplification	X1.0 - X2.0	



Optional Built-in Battery Configuration, with Battery Life Exceeding 4+ Hours.



Brand-Customized UI Interface (MOQ Required)



Lightweight and Portable: Equipped with a Carrying Handle, with a Device Weight of Only 2.5kg



Built-in Patient Recording System and Equipment Usage Count Recording.