

CH

CMOS+LED Technology

Fiber-free design makes it more durable.



High-quality Image

Excellent image gives more details.



Multi-function Buttons

Improve operation flexibility and accuracy.

- ① Zoom in/out
- ② Capture/record
- ③ Freeze/release



Portable

Equipped with image processor supports external monitor.



Specifications

Model	Insertion Tube	Distal End	Working Channel	Working Length	Field of View	Depth of Field	Bending (Up/Down)
CH-1331	3.1mm	3.2mm	1.2mm	450mm	120°	3-50mm	130°/130°
CH-1349	4.8mm	5.0mm	2.0mm	380mm	120°	3-50mm	175°/130°

Image Processor



10.1" Touch screen



1280*800
HD Resolution




HDMI+ AV
output



Rechargeable battery
supporting 4 working-hour



Detachable 64G
SD card



PT. Wickland Medika Imantara
Medical, Aesthetics & Laboratory Supplier

Jl. Palem Sememi Raya No 12 - 12A Kel. Sememi,
Kec. Benowo Surabaya 60198 | 031-33100412