

TH

HD image



TH-E12

Series Video Pleuroscope

Rotation

Improve operation flexibility and accuracy.



Wide bending angle

Solve the problem of limited range of motion between the ribs. The range of observation is much wider.



Angle lock

Bending angle locking function for easy diagnosis and treatment.

(TH-E12 Plus Only)




Suction valve retainer

Unique design prevents suction valve from falling off.



Slim diameter and large working channel

Can match various endoscopic accessories for treatment.



Specifications

Model	Distal End	Insertion Tube	Working Channel	Working Length	Field of View	Depth of Field	Angle Lock	Application
TH-E12 Plus	7.0mm	6.8mm	3.2mm	270mm	120°	3~150mm	Yes	Pneumology
TH-E12	7.0mm	6.8mm	3.2mm	270mm	120°	3~150mm	No	Pneumology

Image Processor



10.1" Touch screen



1280*800 HD resolution



HDMI+AV output



Detachable 64GB SD card



Portable design



Rechargeable battery supporting 4 working-hour



